



See things you've never seen before.

## Dual Mode Fluorescence Module Specifications

**CytoViva™**

### Dimensions

Length = 129mm

Width = 88mm

Height = 144mm

Base = 140mm x 140mm

### Weight

1.75kg ( 62oz.)

### Filter Compatibility

Standard 25mm excitation filters. Emission filter dimensions dependent on microscope requirements.

### Light Guide


Supports 5mm OD (3mm bundle) liquid light guide.

### Light Source

The CytoViva DMF module has been designed to be compatible with the Exfo X-Cite 120 fluorescent light source.

### Warranty

One (1) year from date of purchase.



**Dual Mode Fluorescence Module**  
pictured with CytoViva High Resolution Adapter

