



See things you've never seen before.

Dual Mode Fluorescence Module Specifications

CytoViva™

Dimensions

Length = 129mm

Width = 88mm

Height = 144mm

Base = 140mm x 140mm

Weight

1.75kg (62oz.)

Filter Compatibility

Standard 25mm excitation filters. Emission filter dimensions dependent on microscope requirements.

Light Guide


Supports 5mm OD (3mm bundle) liquid light guide.

Light Source

The CytoViva DMF module has been designed to be compatible with the Exfo X-Cite 120 fluorescent light source.

Warranty

One (1) year from date of purchase.



Dual Mode Fluorescence Module
pictured with CytoViva High Resolution Adapter

